

CCT-RN/Paramedic Treatment Guideline 1201/2201

Medical Assessment and Management Procedures (MAMP)

Page 1 of 2

Follow MAMP Protocol 4201, as applicable, with the following modifications:

A. Patient Assessment.

- 1. Airway and Oxygenation Management:
 - a. If external airway is in place prior to arrival, confirm proper placement and secure device prior to transport.
 - b. If patient's airway is not secure and needs secured, establish patent airway using Airway Management Protocol 4901 and/or Advanced Airway Management (RSI) Guideline 1901/2901.

2. Breathing:

- a. If the patient is in respiratory distress and cannot maintain oxygen saturations >90% by utilizing MAMP Protocol 4201, consider securing airway with endotracheal tube using Airway Management Protocol 4901 and/or Advanced Airway Management (RSI) Guideline 1901/2901.
- b. Assess and document any ventilator settings if applicable.

3. Circulation:

a. Assess vital signs, skin color, temperature, cap refill, vital signs and document.

4. Neuro

- a. Assess level of consciousness, GCS
- b. Assess for any focal neurological deficits.

B. Treatment.

- 1. Initiate or continue at least one (preferably two) IVs of 0.9% normal saline.
- 2. If interfacility transport, assess and document any continuous IV infusions and their rates. Assess to ensure enough meds to complete the transport.
- 3. Monitor ECG, pulse oximeter, end tidal CO2, vital signs at least every 15 minutes, and document.
- 4. If interfacility transport, ensure that any available medical records, x-rays, recent labs are transported with the patient.



CCT-RN/Paramedic Treatment Guideline 1201/2201

Medical Assessment and Management Procedures (MAMP)

Page 2 of 2

- 5. If interfacility transport and prolonged transport time is contemplated, consider urinary catheter insertion.
- 6. If interfacility transport, provide a short update/visit with family at sending facility if time and condition allows, to alleviate anxiety of patient and family.

C. Contact Medical Command enroute with patient report and ETA.

